

PIERTABLE

OR ISLAND BENCH

PIER OR ISLAND TABLES/BENCHES are a great furniture choice for educational science environments. They provide easy access to services and an abundant amount of counter space for each student. Pier tables fit 2-4 students, and are commonly configured in a perimeter layout or centrally clustered and positioned to fit the working requirements of a particular class. Traditional island benches however, are designed to accommodate eight to twelve students, and large enough to provide storage facilities at each lab station.

FEATURES

Solid Wood Construction

Available In 5+ Finishes


Chemical Resistant, UV Finish

Work Surface in Epoxy or Phenolic Resin


Metal Wire Pulls Standard

Optional Center Sink(S)

Drawer, Cabinet or Combination Under Storage



36" H PIER TABLE



34" H ADA PIER TABLE

LAB TABLES

SOLID WOOD FRAME CONSTRUCTION [CHEMICAL RESISTANT - UV FINISH]

5+ WOOD FINISH OPTIONS

FIXED STYLE OR MOBILE BASE

DURABLE WORK SURFACE OPTIONS [EPOXY OR PHENOLIC RESIN]

METAL WIRE PULLS STANDARD

PROTECTIVE RUBBER BOOTS

HEAVY-DUTY LOCKABLE CASTERS



TABLE FRAMES



STEEL COLORS [60% Gloss]

WHITE | BLACK | STORM GRAY | BEIGE | BLUE | DARK BLUE



WORK SURFACES



PHENOLIC RESIN COLORS

TITAN WHITE | BLACK



PLASTIC LAMINATE COLORS

WHITE | BLACK | CLOUDY GRAY | BEIGE | STORM GRAY

DIMENSIONS

ITEM #	TABLETOP	HEIGHT ADJUST	FRAME	WORK SURFACE OPTIONS
LS-KY2424	24"D X 24"L	30" to 36"H	Steel	Phenolic Resin // Plastic Laminate
LS-KY2430	24"D X 30"L	30" to 36"H	Steel	Phenolic Resin // Plastic Laminate
LS-KY2436	24"D X 36"L	30" to 36"H	Steel	Phenolic Resin // Plastic Laminate
LS-KY2448	24"D X 48"L	30" to 36"H	Steel	Phenolic Resin // Plastic Laminate
LS-KY2460	24"D X 60"L	30" to 36"H	Steel	Phenolic Resin // Plastic Laminate
LS-KY2472	24"D X 72"L	30" to 36"H	Steel	Phenolic Resin // Plastic Laminate
LS-KY2496	24"D X 96"L	30" to 36"H	Steel	Phenolic Resin // Plastic Laminate

▶ ALSO AVAILABLE 30" | 36" | 48"D