

MOBILE STORAGE

Mobile Microscope Storage Cabinet

- Features two double locking doors (two per side) that provide secure, safe storage for 24 standard microscopes (6-3/4" x 10" x 15-1/2") but contain removable wire dividers to create a flexible solution for any size microscope
- Constructed of solid oak or maple hardwood and veneers with an earth-friendly UV finish
- Top: 1-1/4" black, high pressure laminate (HPL) with t-mold edge banding
- Rolls easily on heavy-duty swivel casters with locking brakes
- Overall dimensions: 48"W x 24"D x 40-3/4"H

ITEM #
OAK / MAPLE

4701K / 4701M

PRODUCT

Mobile Microscope Storage Cabinet



Large Microscope Storage Cabinet

- Features one large, double locking door cabinet that stores up to fifteen large (9"W x 15"D x 16-1/2") microscopes but contain removable wire dividers to create a flexible solution for any size microscope
- Constructed of solid oak or maple hardwood and veneers with an earth-friendly UV finish
- Rolls easily on heavy-duty swivel casters with locking brakes
- Unit has a finished top, sides and back
- Overall dimensions: 54"W x 24"D x 61-3/4"H

ITEM #
OAK / MAPLE

4711K / 4711M

PRODUCT

Large Microscope Storage Cabinet

Mobile Micro-Charge Station

- Equipped to recharge microscopes safely behind locking doors on each side
- Holds 20 microscopes (7"W x 9"D x 15"H)
- Recharges microscopes with only one electrical connection from the cabinet's side
- Top: 1-1/4" black, high pressure laminate (HPL) with t-mold edge banding
- Constructed of solid oak or maple hardwood and veneers with an earth-friendly UV finish
- Rolls easily on heavy-duty swivel casters with locking brakes
- Overall dimensions: 48"W x 30"D x 42-1/4"H

ITEM #
OAK / MAPLE

4741K / 4741M

PRODUCT

Mobile Microscope Storage Cabinet





Mobile Balance Storage Cabinet

- A safe, efficient, and mobile way to store balances with sliding, locking doors on two sides that provide easy access to the charging station
- Supplied with four adjustable shelves (weight capacity of 40 lbs. per sq. ft.)
- Mobile capacity of 500 lbs. that rolls easily on heavy-duty swivel casters with brakes
- Constructed of solid oak or maple hardwood and veneers with an earth-friendly UV finish
- Top: 1-1/4" black, high pressure laminate (HPL) with t-mold edge banding
- Overall dimensions: 48"W x 24"D x 40"H

ITEM # OAK / MAPLE	PRODUCT
5201K / 5201M	Mobile Balance Storage Cabinet

Tall Mobile Tote Tray Storage

- Models: choose from open-style or locking, two-door style units
- Efficient and mobile tote storage unit that includes 24 classic totes (13-3/4"W x 19"L x 4-3/8"D) for mobile organization
- Choose from oak or maple hardwood and veneers with an earth-friendly UV finish
- Unit has a finished top, back and sides
- Rolls easily on 4" swivel casters
- Overall dimensions: 48-1/2"W x 22-1/2"D x 52"H

ITEM # OAK / MAPLE	PRODUCT
MTTC-4824K / MTTC-4824K	Open-Style Cabinet
MTTC-4824WDK / MTTC-4824WDM	Two-Door Cabinet



Mobile Tote Tray Cabinet

- Unit with a multi-purpose, working surface that includes 21 classic totes (13-3/4"W x 19"L x 4-3/8"D) for mobile organization
- Top: 1-1/4" black, high pressure laminate (HPL) with t-mold edge banding
- Rolls easily on 4" swivel casters
- Choose from oak or maple hardwood and veneers with an earth-friendly UV finish
- Unit has a finished back and sides
- Overall dimensions: 48"W x 24"D x 41-1/2"H

ITEM # OAK / MAPLE	PRODUCT
4751K / 4751M	Mobile Tote Tray Storage

NEW SERIES Duo Mobile Tote Storage



DE-30M3



DE-41K4



DE-65M2

- Choose from three heights: 32-1/2"H, 43-1/2"H, or 66-1/2"H
- Duo-faced units feature totes on two sides
- Rolls easily on 4" swivel casters
- Choose from oak or maple hardwood and veneers with an earth-friendly UV finish
- DE-30 dimensions: 21"W x 37"D x 32-1/2"H with 12 single totes (18-1/2"W x 16-3/4"D x 3"H)
- DE-41 dimensions: 21"W x 37"D x 43-1/2"H with 18 single totes (18-1/2"W x 16-3/4"D x 3"H)
- DE-65 dimensions: 21"W x 37"D x 66-1/2"H with 30 single totes (18-1/2"W x 16-3/4"D x 3"H)
- Unit has a finished top and sides

- Features our Euro-Series (high-impact polystyrene) translucent totes: clear, royal, carrot, or lime
- Clear sturdy glide and tilt tote tray runners secure totes from falling out

STYLE OAK / MAPLE	TOTE OPTION
DE-30K / DE-30M	Euro Translucent Clear (1)
DE-41K / DE-41M	Euro Translucent Royal (2)
DE-65K / DE-65M	Euro Translucent Carrot (3)
	Euro Translucent Lime (4)



Euro Series
Translucent Clear Tote (1)



Euro Series
Royal Translucent Tote (2)



Euro Series
Carrot Translucent Tote (3)



Euro Series
Lime Translucent Tote (4)



NEW SERIES XL Mobile Multi-Tote Storage

- Models: choose from open-style or locking, two-door cabinet units
- Rolls easily on 4" swivel casters
- Choose from oak or maple hardwood and veneers with an earth-friendly UV finish
- Includes nine single totes (12-1/4"W x 16-3/4"L x 3"H); nine double totes (12-1/4"W x 16-3/4"L x 6"H); three quad totes (12-1/4"W x 16-3/4"L x 12"H)
- Features our Euro-Series (high-impact polystyrene) translucent totes: clear, royal, carrot, or lime
- Clear sturdy glide and tilt tote tray runners secure totes from falling out
- Door unit has a three-point locking handle system

- Unit has a finished top, back and sides
- Overall dimensions: 42"W x 21"D x 53-3/4"H

STYLE OAK / MAPLE	TOTE OPTION
(Open) TS-4221K / TS-4221M	Euro Translucent Clear (1)
(Door) TS-4221WDK / TS-4221WDM	Euro Translucent Royal (2)
	Euro Translucent Carrot (3)
	Euro Translucent Lime (4)



NEW SERIES Compact Mobile Tote Storage

- Features locking, door cabinets
- Choose from two heights: 24-3/4"H or 30-3/4"H
- Rolls easily on 3" swivel casters
- Choose from oak or maple hardwood and veneers with an earth-friendly UV finish
- MET-24 unit dimensions are 21"W x 20-1/2"D x 24-3/4"H, with two single totes (18-1/2"W x 16-3/4"L x 3"H), and one double tote (18-1/2"W x 16-3/4"L x 6"H)
- MET-30 unit dimensions are 21"W x 20-1/2"D x 30-3/4"H, with six single totes (18-1/2"W x 16-3/4"L x 3"H)
- Features our Euro-Series (high-impact polystyrene) translucent totes: clear, royal, carrot, or lime

- Clear sturdy glide and tilt tote tray runners secure totes from falling out
- Unit has a finished top, back and sides

STYLE OAK / MAPLE	TOTE OPTION
MET-24K / MET-24M	Euro Translucent Clear (1)
MET-30WDK / MET-30WDM	Euro Translucent Royal (2)
	Euro Translucent Carrot (3)
	Euro Translucent Lime (4)



MTTE-4324M4



MTTE-4324WDK2

NEW SERIES

XL Mobile Single-Tote Storage

- Models: choose from open-style or locking, two-door cabinet units
- Rolls easily on 4" swivel casters
- Choose from oak or maple hardwood and veneers with an earth-friendly UV finish
- Unit has a finished top, back and sides
- Includes 27 single totes (12 1/4" W x 16 3/4" L x 3" H)
- Features our Euro-Series (high-impact polystyrene) translucent totes: clear, royal, carrot or lime
- Clear sturdy glide and tilt tote tray runners secure totes from falling out
- Overall dimensions: 43-1/2"W x 22-1/2"D x 52"H

STYLE OAK / MAPLE	TOTE OPTION
(Open) MTTE-4324K / MTTE-4324M	Euro Translucent Clear (1)
(Door) MTTE-4324WDK / MTTE-4324WDM	Euro Translucent Royal (2)
	Euro Translucent Carrot (3)
	Euro Translucent Lime (4)



TW-4221M4



TW-4221WDK3

NEW SERIES

Mobile Tote & More Storage

- Models: choose from open-style or locking, two-door cabinet units
- Rolls easily on 4" swivel casters
- Choose from oak or maple hardwood and veneers with an earth-friendly UV finish
- Unit has a finished top, back and sides
- Includes six single totes (12-1/4"W x 16-3/4"L x 3"H), six double totes (12-1/4"W x 16-3/4"L x 6"H), and two quad totes (12-1/4"W x 16-3/4"L x 12"H)
- Includes backpack/coat space (13"W x 37"H) with interior hook
- Features our Euro-Series (high-impact polystyrene) translucent totes: clear, royal, carrot, or lime
- Clear sturdy glide and tilt tote tray runners secure totes from falling out

- Door unit has a three-point locking handle system
- Overall dimensions: 42"W x 21"D x 53-3/4"H

STYLE OAK / MAPLE	TOTE OPTION
(Open) TW-4221K / TW-4221M	Euro Translucent Clear (1)
(Door) TW-4221WDK / TW-4221WDM	Euro Translucent Royal (2)
	Euro Translucent Carrot (3)
	Euro Translucent Lime (4)

NEW OPTIONS

Open & Closed Tote Cabinets



- Models: choose from open-style or locking, two-door cabinet units
- Rolls easily on 4" swivel casters
- Choose from oak or maple hardwood and veneers with an earth-friendly UV finish
- Unit has a finished back and sides
- Top: 1-1/4" black, high pressure laminate (HPL) with t-mold edge banding
- 4801 dimensions are 48"W x 24"D x 40-1/2"H; with 12 single totes (12-1/4"W x 16-3/4"D x 3"H); 3 double totes (12-1/4"W x 16-3/4"D x 6"H); and one large wood drawer (39"W x 20"D x 4"H)
- 4811 dimensions are: 48"W x 24"D x 43-1/2"H; with 27 single totes (12-1/4"W x 16-3/4"D x 3"H)
- Features Euro-Series (high-impact polystyrene) translucent totes: clear, royal, carrot, or lime
- Clear sturdy glide and tilt tote tray runners that secure totes from falling out

STYLE OAK / MAPLE	TOTE OPTION
(Open) 4801K / 4801M	Euro Translucent Clear (1)
(Door) 4811K / 4811M	Euro Translucent Royal (2)
	Euro Translucent Carrot (3)
	Euro Translucent Lime (4)



Clear
Translucent
Totes



Royal
Translucent
Totes



Carrot
Translucent
Totes



Lime
Translucent
Totes

NEW OPTIONS

Write-n-Roll Storage



- Models: choose from Classic or Euro-Series tote series units
- Rolls easily on 4" swivel casters
- Choose from oak or maple hardwood and veneers with an earth-friendly UV finish
- Adjustable, 3/4" white markerboard top with a 3mm black PVC edge,
- Features two 10-position arms for smooth adjustability with soft-close "no pinch" safety feature
- Additional storage under the top
- 4901 dimensions: (36"W x 24"D x 35-1/2"H) with five single (classic) totes (13-3/4"W x 19"L x 4-3/8"D); and three wood drawers (13-1/2"W x 16-3/4"D x 2-1/8"H)
- 4821 dimensions: (36"W x 24"D x 35-1/2"H) with six single Euro totes (12-1/4"W x 16-3/4"D x 3"H); and three wood drawers (14-5/8"W x 16-3/4"D x 2-1/8"H)
- 4821 series features our Euro-Series (high-impact polystyrene) translucent totes: clear, royal, carrot, or lime; with clear sturdy glide and tilt tote tray runners that secure totes from falling out

STYLE OAK / MAPLE	TOTE OPTION
4901K / 4901M	Classic Tan

STYLE OAK / MAPLE	TOTE OPTION
4821K / 4821M	Euro Translucent Clear (1)
	Euro Translucent Royal (2)
	Euro Translucent Carrot (3)
	Euro Translucent Lime (4)

Economy Mobile Storage Cabinet

- Model: locking, two-door cabinet unit
- Choose from oak or maple hardwood and veneers with an earth-friendly UV finish
- Top: 1-1/4" black, high pressure laminate (HPL) with t-mold edge banding
- Features sliding, locking doors for easy access and secure storage
- Supplied with one, fixed shelf (weight capacity of 40 lbs. per sq. ft.)
- Mobile capacity of 500 lbs. that rolls on heavy-duty swivel casters
- Overall dimensions: 48"W x 24"D x 36"H

ITEM # OAK / MAPLE	TOP
4101K / 4101M	1-1/4"HPL
4102K / 4102M	1-1/4" ChemGuard



Mobile Storage Cabinet

- Model: locking, two-door cabinet unit
- Choose from oak or maple hardwood and veneers with an earth-friendly UV finish
- Top options: 1-1/4" black, high pressure laminate (HPL) or 1-1/4" ChemGuard
- Mobile capacity of 500 lbs. that rolls on heavy-duty swivel casters
- Overall dimensions: 36"W x 24"D x 36"H

ITEM # OAK / MAPLE	TOP
4401K / 4410M	1-1/4"HPL
4402K / 4402M	1-1/4" ChemGuard

Mobile Demo Cart

- Model: open-style cart
- Choose from oak or maple hardwood and veneers with an earth-friendly UV finish
- Top options: 1-1/4" black, high pressure laminate (HPL), 1-1/4" ChemGuard, or 1-3/4" maple butcher block
- Legs are 2-1/4" square and are fastened to the 4-1/2" high apron with factory installed 5/16" hanger bolts through the heavy-duty 13-gauge steel corner brace
- Includes leg braces with an auxiliary, lower plywood shelf
- Mobile capacity of 500 lbs. that rolls on heavy-duty swivel casters
- Overall dimensions: 36"W x 24"D x 36"H

ITEM # OAK / MAPLE	TOP
4501K / 4501M	1-1/4"HPL
4502K / 4502M	1-1/4" ChemGuard
PDC-2436K / PDC-2436M	Maple BB

