



## TALL STORAGE

Tall storage cabinets are 82" high and depths of 16" and 22". Units have four adjustable shelves and center shelf fixed. All dimensions are nominal.



## Open Storage




### 16" DEEP

T188216-000 - 18"W  
T248216-000 - 24"W

### 22" DEEP

T188222-000 - 18"W  
T248222-000 - 24"W





### 16" DEEP

T308216-000 - 30"W  
T358216-000 - 35"W  
T428216-000 - 42"W  
T478216-000 - 47"W

### 22" DEEP

T308222-000 - 30"W  
T358222-000 - 35"W  
T428222-000 - 42"W  
T478222-000 - 47"W

## Solid Wood Door




### 16" DEEP

T188216-100 - 18"W RH  
T188216-101 - 18"W LH  
T248216-100 - 24"W RH  
T248216-101 - 24"W LH

### 22" DEEP

T188222-100 - 18"W RH  
T188222-101 - 18"W LH  
T248222-100 - 24"W RH  
T248222-101 - 24"W LH





### 16" DEEP

T308216-200 - 30"W  
T358216-200 - 35"W  
T428216-200 - 42"W  
T478216-200 - 47"W

### 22" DEEP

T308222-200 - 30"W  
T358222-200 - 35"W  
T428222-200 - 42"W  
T478222-200 - 47"W

## Framed Glass Door



### 16" DEEP

T188216-150 - 18"W RH  
T188216-151 - 18"W LH  
T248216-150 - 24"W RH  
T248216-151 - 24"W LH

### 22" DEEP

T188222-150 - 18"W RH  
T188222-151 - 18"W LH  
T248222-150 - 24"W RH  
T248222-151 - 24"W LH




### 16" DEEP

T308216-250 - 30"W  
T358216-250 - 35"W  
T428216-250 - 42"W  
T478216-250 - 47"W


### 22" DEEP

T308222-250 - 30"W  
T358222-250 - 35"W  
T428222-250 - 42"W  
T478222-250 - 47"W

Special use tall cases are 82" high and 16" deep. All dimensions are nominal.




## Wardrobe Cabinets



**T248216-120 - 24"W RH**  
**T248216-121 - 24"W LH**

This unit features a garment hangar with pull-out carrier, hat shelf, mirror and solid door.


## Notebook Cabinets



**T248216-130 - 24"W RH**  
**T248216-131 - 24"W LH**

Horizontal compartment for storing flat materials including notebooks. (29) individual compartments each 2 3/8" high.


## Skeleton Cabinets



**T248216-140 - 24"W RH**  
**T248216-141 - 24"W LH**

Comes with framed glass door and pull-out carrier for skeleton.


## Torso Cabinets



**T248216-160 - 24"W RH**  
**T248216-161 - 24"W LH**


Comes with framed glass door and two adjustable shelves and pull-out shelf for top.

## Microscope Cabinets



**T248216-170 - 24"W RH**  
**T248216-171 - 24"W LH**


Accommodates up to (15) microscopes.



**T358216-270 - 35"W**  
**T478216-270 - 47"W**

Unit stores as many as (30) microscopes. Hinged glass doors allow for easy access, as well as instant recognition.


## Tubing Cabinets



**T248216-180 - 24"W RH**  
**T248216-181 - 24"W LH**

Slanted storage compartments allow easy access to glass tubing.


## Axis 3<sup>®</sup> Storage Cabinets




**T478216-799 - 47"W**

Unit features framed glass doors with lock. Used for storage of wireless keyboards and various accessories for the Axis 3<sup>®</sup> student table.

Special use tall cases on this page are 82" high and 22" deep.  
All dimensions are nominal.




## Display



**T358222-700 - 35"W**  
**T478222-700 - 47"W**

Provided with (4) adjustable shelves, one fixed shelf and framed glass doors.


## Compartmental



**T358222-712 - 35"W**  
**T478222-712 - 47"W**

Three vertical spaces and nine movable shelves. Permits arrangement of compartments front to back.


## Blind Corner



**T478222-100 - 47"W RH**  
**T478222-101 - 47"W LH**


Two full length adjustable shelves and one fixed shelf. Hinged doors. RH version shown.

## Instructor




**T478222-710 - 47"W RH**  
**T478222-711 - 47"W LH**

Cabinet includes (2) file drawers, (3) sloping display shelves, mirror and a wardrobe compartment with clothes rod & fixed hat shelf. RH version shown.



**T358222-720 - 35"W**  
**T478222-720 - 47"W**


Left hand section fitted with (5) adjustable shelves. Right hand section with fixed shelf, clothes rod & mirror on door.



**T358222-721 - 35"W**  
**T478222-721 - 47"W**

Right hand section fitted with (5) adjustable shelves. Left hand section with fixed shelf, clothes rod & mirror door.

## Tote Tray - 36




**T478222-236 - 47"W**

Unit comes complete with (36) tote trays.

**Tote Tray Size:**  
19 1/2"L X 13 1/4" W X 4 1/2"H - O.D.

## Tote Tray - 50




**T478222-210 - 47"W**

Unit comes complete with (50) tote trays.

**Tote Tray Size:**  
19 1/2"L X 8 13/16" W X 4 1/2"H - O.D.

## Tote Tray - 60




**T478222-220 - 47"W**

Unit comes complete with (60) tote trays.

**Tote Tray Size:**  
19 1/2"L X 8 13/16" W X 4 1/2"H - O.D.


## Carryall



**T478222-730 - 47"W**

The Sheldon Carryall & Carryall Storage make possible the most economical & easily controlled student equipment storage. This unit includes (28) Carryalls.


## Apparatus Panel Storage



**T448222-750 - 44"W**

Upper section holds up to (28) blank panel boards & approximately (14) ready-mounted apparatus panel boards in the vertical position. Lower section with (3) tray drawers to accommodate loose apparatus & holds up to (16) blank panel boards & approximately (10) ready mounted apparatus panel boards in the horizontal position.

## Tool



**T478222-740 - 47"W**

Pegboard lining upper area for hooks to store tools. Lower portion of cabinet with two adjustable shelves.